

Manawa

URETERAL ACCESS SHEATH

SMOOTH INSERTION
FOR CONTINUOUS
URETERAL ACCESS




Multiple configuration available

The Manawa ureteral access sheath provides ureteral dilation with continuous working channel for the introduction of ureteroscopes and instruments during ureteroscopic procedure.

Effective access for safer procedure

- Protect the ureter during repeated instrument exchange minimizing the risk of trauma
- The PTFE lined sheath protects delicate instruments and ureteroscopes from damage
- Faster procedure time and enhanced wash-out of stone fragments
- Reduce the need of balloon dilators

Smooth insertion and gentle dilation



Order Number	Internal Diameter French Size	External Diameter French Size	Length
<i>Manawa access sheath</i>			
UAS-0109520S	9.5 Fr	11.5 Fr	20 cm
UAS-0109528S	9.5 Fr	11.5 Fr	28 cm
UAS-0109535S	9.5 Fr	11.5 Fr	35 cm
UAS-0109545S	9.5 Fr	11.5 Fr	45 cm
UAS-0110720S	10.7 Fr	12.7 Fr	20cm
UAS-0110728S	10.7 Fr	12.7 Fr	28cm
UAS-0110735S	10.7 Fr	12.7 Fr	35cm
UAS-0110745S	10.7 Fr	12.7 Fr	45cm
UAS-0112020S	12.0 Fr	14.0 Fr	20 cm
UAS-0112028S	12.0 Fr	14.0 Fr	28 cm
UAS-0112035S	12.0 Fr	14.0 Fr	35 cm
UAS-0112045S	12.0 Fr	14.0 Fr	45 cm
UAS-0114020S	14.0 Fr	16.0 Fr	20 cm
UAS-0114028S	14.0 Fr	16.0 Fr	28 cm
UAS-0114035S	14.0 Fr	16.0 Fr	35 cm
UAS-0114045S	14.0 Fr	16.0 Fr	45 cm